

The Intelsius BIOSHIELD range of unique patented systems designed to meet handling processes from sample extraction to laboratory presentation. BIOSHIELDS are a highly protective single-use shipping solution which are available for a variety of sample capacities, helping reduce waste and cut costs associated with sample transportation



The systems are well suited to a wide variety of samples, including solids, liquids, and samples held in fixative. The BIOSHIELD also ensures that all samples will arrive at the same time for a simple and efficient sample reception, helping to reduce total costs.

Example Components Supplied:

- Outer Carton printed in compliance with transport regulations
- 100ml Pots, 95kPa IATA compliant
- Pot Tray
- Leak-proof 95kPa secondary container
- Security Seal Label
- Techni-Ice Gel Sheet (Next Day Only)

BIOSHIELD 'Same Day' RANGE

PRODUCT	PRODUCT CODE (INC. POTS)	PRODUCT CODE (EXC. POTS)	PRODUCT CODE (INC. POTS, SPOONS & FORCEPS)	SAMPLE CAPACITY	EXTERNAL DIMENSIONS	VOLUMETRIC WEIGHT
BIOSHIELD 6	BSE061	TBC	BSE055	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BIOSHIELD 10	TBC	TBC	BSE038	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BIOSHIELD SR 20	BSE056	BSE065	BSE057	20	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BIOSHIELD SR 40	BSE058	BSE066	BSE059	40	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BIOSHIELD TR 80	TBC	BSE021	TBC	80	495 x 355 x 245 mm 19.5 x 14 x 9.7 in	7.2 kgvol ^c 15.9 lbsvol ^c

BIOSHIELD 'Next Day' RANGE

PRODUCT	PRODUCT CODE (INC. POTS)	PRODUCT CODE (EXC. POTS)	PRODUCT CODE (INC. POTS, SPOONS & FORCEPS)	SAMPLE CAPACITY	EXTERNAL DIMENSIONS	VOLUMETRIC WEIGHT
BIOSHIELD 6	BSE062	TBC	BSE051	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BIOSHIELD 10	TBC	TBC	TBC	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BIOSHIELD SR 20	BSE050	BSE063	BSE064	20	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BIOSHIELD SR 40	BSE041	BSE042	BSE067	40	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BIOSHIELD TR 80	BSE019	BSE007	TBC	80	495 x 355 x 245 mm 19.5 x 14 x 9.7 in	7.2 kgvol ^c 15.9 lbsvol ^c