

The BIOTHERM Category A 2-8°C shippers ensure that consignments are fully compliant to IATA, ADR & 49 CFR (DOT) transport regulations. The systems are pre-qualified for 96 hours against the demanding Intelcius High Performance Summer and Winter Temperature Profiles.



The systems are designed to optimise shipping volumes and minimise freight costs. Each features unique sturdy handles to enable precise handling throughout the cold chain, and a tongue and groove lid and base to enable quick and simple assembly. Pack-out instruction sheets are supplied to ensure correct assembly thereby ensuring peak performance, and a security seal label ensures tamper evidence and security confidence.

#### Components Supplied:

- Pack-out Instruction Sheet
- 2 Litre PATHOPAK 95kPa IATA compliant leak-proof secondary container
- Payload Carton
- Absorbent Material
- Qualified 2-8°C Kit
- Molded Carrier Handle
- EPS Insulated Shipper
- Label Pack including: Infectious Class Label, Dry Ice, Overpack, Security Seal Label

PRODUCT	CODE	SECONDARY CONTAINER INTERNAL DIAMETER x HEIGHT	EXTERNAL DIMENSIONS	VOLUMETRIC WEIGHT	MASS	DURATION (HOURS)
BIOTHERM 30A 2-8° HP 96 S	BT056	125 x 170 mm 4.9 x 6.6 in	465 x 458 x 478 mm 18.3 x 18.1 x 18.9 in	17 kgvol <sup>c</sup> 37.4 lbsvol <sup>c</sup>	10.2 kg 22.4 lbs	96
BIOTHERM 30A 2-8° HP 96 W	BT057	125 x 170 mm 4.9 x 6.6 in	465 x 458 x 478 mm 18.3 x 18.1 x 18.9 in	17 kgvol <sup>c</sup> 37.4 lbsvol <sup>c</sup>	9 kg 19.8 lbs	96

**Product Name Key:** [Range] [EPS Size & Shipping Category] [Temperature] [Test Profile] [Duration] [Season]