

The PATHOSHIELD 4 is UN Certified for Category A shipments, meeting IATA (Air), ADR (Euro Road) and 49 CFR transport regulations. The system is designed for shipping a wide range of primaries, including blood vials, swabs and specimen containers.



The PATHOSHIELD 4 Category A Shipper is designed to optimise shipping volume and minimise freight costs. It can be flat-packed for efficient storage or shipping prior to use. The security seal label provides tamper evidence and security confidence, whilst the A4 PATHOPOUCH 95kPa flexible secondary container ensures compliance to the required shipping regulations.

#### Components Supplied:

- PATHOSHIELD P4 Outer Carton
- Flexible 95kPa compliant secondary container. (A4 PATHOPOUCH Bag)
- 6 Bay Absorbent Material
- Security Seal, and Infectious Specimen Labels

PRODUCT	CASE QTY	CODE	POUCH OPENING SIZE	POUCH INTERNAL DIMENSIONS	EXTERNAL DIMENSIONS	VOLUMETRIC WEIGHT	MASS
PATHOSHIELD 4	25	PA002	220mm / 8.6in	220 x 240 mm 8.6 x 9.4 in	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol <sup>c</sup> 2.1 lbsvol <sup>c</sup>	0.2 kg 0.5 lbs